Equipment Training Sheet - MTS Insight 10

Name:	Position:
Department:	Phone:
e-mail:	Director:
Staff Trainee Initials Initials	
	 Anyone needing access to MSE equipment needs to contact MSE tech staff: For initial training For subsequent usage (will receive equipment card) No secondary training of the equipment
	 Review of hazards, controls, personal protective equipment: Always wear safety glasses when operating the MTS, also be aware of anyone close by to either move to a safe distance or to wear safety glasses Always use machine shield when testing samples in the MTS Pinch and crush hazard due to moving crosshead, extra caution should be taken Be mindful of the location of the emergency stop button
	 Review of standard operating procedure: http://www.mse.engin.umich.edu/internal/procedures/mts-insight-10 Before testing calculate the maximum load conditions of your sample to determine if the machine\load cell is appropriate Make sure the correct load cell is attached to the testing frame before turning it on, 1kN and 10kN load cells available A 2 inch gauge section Instron extensometer is available for use, attach this before testing as well Be aware that certain grips and testing fixtures have lower limits than the load cell: There is a set of fine teethed tensile grips with a limit of 5kN There is a set of 3&4 point bending testing fixtures with a limit of 5kN Remember to set the hardware limits on the testing frame to prevent crashing the grips and causing damage to the crosshead and\or load cell
	Available supplies to help in the procedure: • Rulers and calipers • Tweezers • Safety glasses • Several different pairs of grips/fixtures • 2" gauge extensometer
Staff Signature/da	ate:

Trainee Signature/date: