

Equipment Training Sheet - Bede D1

Name: _____

Position: _____

Department: _____

Phone: _____

e-mail: _____

Director: _____

Staff Initials	Trainee Initials	
		<p><u>Anyone needing access to MSE equipment needs to contact MSE tech staff:</u></p> <ul style="list-style-type: none"> • For initial training • For subsequent usage (will receive equipment card) • No secondary training of the equipment
		<p><u>Review of hazards, controls, personal protective equipment:</u></p> <ul style="list-style-type: none"> • Machine generates X-rays, Take all necessary precautions to avoid exposure to x-rays. • Wear dosimeter rings at all times while operating Bede D1 to monitor any radiation exposure; Store dosimeter in the cabinet after you are done with Bede. • Bring the MSDS for the material to be studied • This machine uses high voltages, do not disassemble or touch inside components.
		<p><u>Review of standard operating procedure:</u></p> <ul style="list-style-type: none"> • http://www.mse.engin.umich.edu/internal/procedures/bede-d1 • Report any problem with the equipment to Staff. If you think the equipment has failed and do not attempt to repair it yourself. • Usage time begins when you occupy this instrument.
		<p><u>Available supplies to help in the procedure:</u></p> <ul style="list-style-type: none"> • Plastic Tweezers, gloves, Scissors • Double sided tape

Staff Signature/date: _____

Trainee Signature/date: _____